

# Rapid Lambda Control Experiment (SQK-RAD004)

Version: RSE\_9045\_v1\_revY\_14Aug2019  
 Last update: 22/12/2021



Flow Cell Number: .....

DNA Samples: .....

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> Rapid Sequencing Kit (SQK-RAD004)	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Microfuge
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP002)	<input type="checkbox"/> 0.2 ml thin-walled PCR tubes	<input type="checkbox"/> Timer
	<input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, cat # AM9937)	<input type="checkbox"/> Thermal cycler or heat block at 30°C and 80°C
		<input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P1000

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><b>Library preparation</b></p> <p>DNA tagmentation</p> <p>Thaw kit components at RT, spin down briefly using a microfuge and mix by pipetting as indicated by the table below:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Lambda DNA (50 µg/ml) (LMD): thaw at RT, briefly spin down, mix well by pipetting</li> <li><input type="checkbox"/> Fragmentation Mix (FRA): not frozen, briefly spin down, mix well by pipetting</li> <li><input type="checkbox"/> Rapid Adapter (RAP): not frozen, briefly spin down, mix well by pipetting</li> <li><input type="checkbox"/> Sequencing Buffer (SQB): thaw at RT, briefly spin down, mix well by pipetting*</li> <li><input type="checkbox"/> Loading Beads (LB): thaw at RT, briefly spin down, mix by pipetting or vortexing immediately before use</li> <li><input type="checkbox"/> Flush Buffer (FLB): thaw at RT, briefly spin down, mix well by pipetting*</li> <li><input type="checkbox"/> Flush Tether (FLT): thaw at RT, briefly spin down, mix well by pipetting</li> </ul> <p><input type="checkbox"/> Once thawed, keep all the kit components on ice.</p> <p>In a 0.2 ml thin-walled PCR tube, mix the following:</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> 7.5 µl Lambda DNA</li> <li><input type="checkbox"/> 2.5 µl FRA</li> </ul> <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Incubate the tube at 30°C for 1 minute and then at 80°C for 1 minute. Briefly put the tube on ice to cool it down.</p>	
<p>375 ng tagmented Lambda DNA in 10 µl is taken into the next step.</p>	
<p><b>Adapter attachment</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Add 1 µl of RAP to the tube.</li> <li><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</li> </ul>	

# Rapid Lambda Control Experiment (SQK-RAD004)

Version: RSE\_9045\_v1\_revY\_14Aug2019  
 Last update: 22/12/2021



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Incubate the reaction for 5 minutes at RT.	
The prepared DNA library is used for loading into the flow cell. Store the library on ice until ready to load.	
<b>Priming and loading the SpotON Flow Cell</b>	
<p><b>IMPORTANT</b></p> <input type="checkbox"/> Please note that the Sequencing Tether (SQT) tube will NOT be used in this protocol. It is provided in the kit for potential future product compatibility.	
<input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of Flush Buffer (FB) at RT.  <input type="checkbox"/> Mix the Sequencing Buffer (SQB), Flush Tether (FLT) and Flush Buffer (FB) tubes by vortexing and spin down at RT.  <input type="checkbox"/> Open the MinION Mk1B lid and slide the flow cell under the clip.  <input type="checkbox"/> Slide the priming port cover clockwise to open the priming port.	
<p><b>IMPORTANT</b></p> <input type="checkbox"/> Take care when drawing back buffer from the flow cell. Do not remove more than 20-30 $\mu$ l, and make sure that the array of pores are covered by buffer at all times. Introducing air bubbles into the array can irreversibly damage pores.	
<p>After opening the priming port, check for a small air bubble under the cover. Draw back a small volume to remove any bubbles (a few <math>\mu</math>l):</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Set a P1000 pipette to 200 <math>\mu</math>l</li> <li><input type="checkbox"/> Insert the tip into the priming port</li> <li><input type="checkbox"/> Turn the wheel until the dial shows 220-230 <math>\mu</math>l, or until you can see a small volume of buffer entering the pipette tip</li> </ul> <p>Note: Visually check that there is continuous buffer from the priming port across the sensor array.</p> <input type="checkbox"/> To prepare the flow cell priming mix, add 30 $\mu$ l of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FB), and mix by vortexing at RT.  <input type="checkbox"/> Load 800 $\mu$ l of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes. During this time, prepare the library for loading by following the steps below.  <input type="checkbox"/> Thoroughly mix the contents of the Loading Beads (LB) tubes by vortexing.  In a new tube, prepare the library for loading as follows: <ul style="list-style-type: none"> <li><input type="checkbox"/> 34 <math>\mu</math>l Sequencing Buffer (SQB)</li> <li><input type="checkbox"/> 25.5 <math>\mu</math>l Loading Beads (LB), mixed immediately before use</li> <li><input type="checkbox"/> 4.5 <math>\mu</math>l Nuclease-free water</li> <li><input type="checkbox"/> 11 <math>\mu</math>l DNA library</li> </ul>	
<p><b>IMPORTANT</b></p> <input type="checkbox"/> The Loading Beads (LB) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use.	

# Rapid Lambda Control Experiment (SQK-RAD004)

Version: RSE\_9045\_v1\_revY\_14Aug2019  
Last update: 22/12/2021



Flow Cell Number: .....

DNA Samples: .....

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Complete the flow cell priming:</p> <ul style="list-style-type: none"><li><input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible.</li><li><input type="checkbox"/> Load 200 µl of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles.</li> <li><input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading.</li> <li><input type="checkbox"/> Add 75 µl of sample to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next.</li> <li><input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION Mk1B lid.</li></ul>	
<p><b>Ending the experiment</b></p>	
<ul style="list-style-type: none"><li><input type="checkbox"/> After your sequencing experiment is complete, if you would like to reuse the flow cell, please follow the Wash Kit instructions and store the washed flow cell at 2-8°C, OR</li> <li><input type="checkbox"/> Follow the returns procedure by washing out the flow cell ready to send back to Oxford Nanopore.</li></ul>	
<p><b>IMPORTANT</b></p> <ul style="list-style-type: none"><li><input type="checkbox"/> If you encounter issues or have questions about your sequencing experiment, please refer to the Troubleshooting Guide that can be found in the online version of this protocol.</li></ul>	